


VERTICAL LAMINAR FLOW CABINET CLB -201-04



Downflow Velocity 0.3m/s ~ 0.4m/s (Variable Set Speed)
Filter Efficiency HEPA Filter ≥ 99.99% at 0.3μm / Class 100
Noise Level ≤ 65dB(A) at Normal Speed
UV Lamp Germicidal UV-Lamp / Auto OFF
Door 5mm Glass Sliding Type
Work Top Stainless Steel SUS304


Air Flow



HORIZONTAL LAMINAR FLOW CABINET CLB-202-04



Air Flow



DUCTLESS STORAGE CABINET



- The casters are attached to install easily.
- There is no need to construct duct so you can move it freely and keep reagent safe as well.
- Alarm let users know the abnormality, when blower would be out of order.
- Spill tray is set under the body that minimizes to out flow harmful reagent.
- The brown colour door is made of transparent quality to protects light from outside and also make users check reagent easily inside of the cabinets.
- Filter test hole is set on the top; you can check the saturation of filtering using gas detection kit.

SBWS (SAFETY TRIMMING WORK STATION) CLE -081F-03



Safety trimming working station can make you easy to wash and protect users from harmful gas, bad smell in animal trimming or test.

Dimensions(External)	1600x525x2160
WxDxH (mm) (Internal)	1520x455x1615
Filter	Carbon Filter (General Type)
Blower	Centrifugal Blower
Noise Level	≤ 65dB (A)
Volume of Air Treated	150m ³ /h
Shelves	8ea
Door	Acrylic 5mm(Brown Colour)

FUME HOOD CLE-101-05



External Frame Epoxy powder coated aluminum
Internal Plate Phenolic laminate
Work top Ceramite
Window Tempered safety glass
 Ergonomic air foil
Receptacle 220V/110V(Selectable)
Cup sink PP
Cabinet Reagent storage cabinet
Switch Lamp/Blower ON/OFF
Cabinet ventilation pipe PVC

FUME HOOD TEST

- Flow Visualization
- Face Velocity Measurement
- Tracer Gas Containment Test

FUME HOOD WALK-IN TYPE CLE-110

